

REF **SAC-01218-B** (IPNXXXXXX)
**18** Ga. **12** cm catheter length **.025** inch dia. spring-wire guide

**ARROW**  
**ErgoPack**®


# Arterial Catheterization

## Contents:

- |  |   |
|--|---|
| 1: Catheter: 18 Ga. x 12 cm Radiopaque PEBA  | 1: Surgical Apparel: Impervious Gown  |
| 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: Dressing: I.V. Advanced 10 cm x 12 cm                                      |
| 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW   | 1: Dressing: BIOPATCH® <sup>3</sup>   |
| 1: Injection Needle: Eclipse™ <sup>1</sup> 23 Ga. x 1-1/2" (3.81 cm)   | 1: Transducer Cover: 14 cm x 147 cm   |
| 1: Injection Needle: Eclipse™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)   | 1: Surgical Apparel: Bouffant Cap   |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)  | 1: Surgical Apparel: Mask with Eye Shield                                     |
| 1: Syringe: 3 mL Luer-Lock   | 1: Suture: 3-0 Silk with Curved Needle  |
| 1: Syringe: 5 mL Luer-Slip   | 1: HemoHopper® Fluid Receptacle   |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep   | <sup>1</sup> A trademark of Becton, Dickinson and Company.                    |
| 2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® <sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint                         | <sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries. |
| 1: Dust Cap: Non-Vented  | <sup>3</sup> A registered trademark of Johnson & Johnson Corporation.         |
| 1: SharpsAway® Disposal Cup  |   |
| 1: SharpsAway® II Locking Disposal Cup   |   |
| 1: Maximal Barrier Drape™ with 4" fenestration   |   |
| 1: Towel   |   |
| 1: Needle Holder   |   |
| 1: Safety Scalpel: #11   |   |
| 1: Checklist Sheet   |   |
| 1: Medication Label: 1% Lidocaine  |   |
| 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)  |   |
| 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)  |   |

**Catheter Priming Volume: 0.21 mL Priming volume is approximate.**

### Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

