

REF **SAC-01618-C** (IPNXXXXXX)
18 Ga. **16** cm catheter length **.025** inch dia. spring-wire guide

ARROW
ErgoPack®


Arterial Catheterization

Contents:

- | | |
|--|---|
| <ul style="list-style-type: none"> 1: Catheter: 18 Ga. x 16 cm Radiopaque PEBA 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW 1: Injection Needle: Eclipse™¹ 23 Ga. x 1-1/2" (3.81 cm) 1: Injection Needle: Eclipse™¹ 25 Ga. x 1" (2.54 cm) 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) 1: Syringe: 3 mL Luer-Lock 1: Syringe: 5 mL Luer-Slip 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep 2: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®² One-Step Solution with Hi-Lite Orange™ Tint 1: Dust Cap: Non-Vented 1: SharpsAway® Disposal Cup 1: SharpsAway® II Locking Disposal Cup 1: Maximal Barrier Drape™ with 4" fenestration 1: Towel 1: Needle Holder 1: Safety Scalpel: #11 1: Checklist Sheet 1: Medication Label: 1% Lidocaine 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) | <ul style="list-style-type: none"> 1: Surgical Apparel: Impervious Gown 1: Dressing: Tegaderm®³ CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm) 1: Transducer Cover: 14 cm x 147 cm 1: Surgical Apparel: Bouffant Cap 1: Surgical Apparel: Mask with Eye Shield 1: Suture: 3-0 Silk with Curved Needle 1: HemoHopper® Fluid Receptacle <p>¹A trademark of Becton, Dickinson and Company.
 ²A registered trademark of CareFusion or one of its subsidiaries.
 ³A registered trademark of 3M Company.</p> |
|--|---|

Catheter Priming Volume: 0.24 mL Priming volume is approximate.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

