

REF **SAC-01218-C** (IPNXXXXXX)
**18** Ga. **12** cm catheter length **.025** inch dia. spring-wire guide

**ARROW**  
**ErgoPack**®


# Arterial Catheterization

## Contents:

- |  |   |
|--|---|
| <ul style="list-style-type: none"> <li>1: Catheter: 18 Ga. x 12 cm Radiopaque PEBA</li> <li>1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer</li> <li>1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW</li> <li>1: Injection Needle: Eclipse™<sup>1</sup> 23 Ga. x 1-1/2" (3.81 cm)</li> <li>1: Injection Needle: Eclipse™<sup>1</sup> 25 Ga. x 1" (2.54 cm)</li> <li>1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)</li> <li>1: Syringe: 3 mL Luer-Lock</li> <li>1: Syringe: 5 mL Luer-Slip</li> <li>1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep</li> <li>2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint</li> <li>1: Dust Cap: Non-Vented</li> <li>1: SharpsAway® Disposal Cup</li> <li>1: SharpsAway® II Locking Disposal Cup</li> <li>1: Maximal Barrier Drape™ with 4" fenestration</li> <li>1: Towel</li> <li>1: Needle Holder</li> <li>1: Safety Scalpel: #11</li> <li>1: Checklist Sheet</li> <li>1: Medication Label: 1% Lidocaine</li> <li>4: Gauze Pad: 2" x 2" (5 cm x 5 cm)</li> <li>4: Gauze Pad: 4" x 4" (10 cm x 10 cm)</li> </ul> | <ul style="list-style-type: none"> <li>1: Surgical Apparel: Impervious Gown</li> <li>1: Dressing: Tegaderm®<sup>3</sup> CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)</li> <li>1: Transducer Cover: 14 cm x 147 cm</li> <li>1: Surgical Apparel: Bouffant Cap</li> <li>1: Surgical Apparel: Mask with Eye Shield</li> <li>1: Suture: 3-0 Silk with Curved Needle</li> <li>1: HemoHopper® Fluid Receptacle</li> </ul> <p><sup>1</sup>A trademark of Becton, Dickinson and Company.<br/> <sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.<br/> <sup>3</sup>A registered trademark of 3M Company.</p> |
|--|---|

---

**Catheter Priming Volume: 0.21 mL Priming volume is approximate.**


---

### Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

