



REF **SAC-01618-1A** (IPNXXXXXX)

18^{Ga.} **16**^{cm} catheter length **.025**^{inch dia.} spring-wire guide

Arterial Catheterization



Contents:

- 1: Catheter: 18 Ga. x 16 cm Radiopaque PEBA Origin: XX XXXX XX XXX
XXXXXXXXXXXXXXXXXXXX

1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer Origin: XX XXXX XX XXX
XXXXXXXXXXXXXXXXXXXX

1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW Origin: XX XXXX XX XXX
XXXXXXXXXXXXXXXXXXXX

1: Injection Needle: Eclipse™ 1 25 Ga. x 1" (2.54 cm) Origin: XX XXXX XX XXX
XXXXXXXXXXXXXXXXXXXX

1: Injection Needle: SafetyGlide™ 1 23 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX
XX XXX XXXXXXXX XXXXXXXXXX

1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX
XXXXXXXXXXXXXXXXXXXX

1: Syringe: 3 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX

1: Syringe: 5 mL Luer-Slip Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX

1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep Origin: XX XXXX XX XXX
XXXXXXXXXXXXXXXXXXXX

2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® 2 One-Step Solution with
Hi-Lite Orange™ Tint Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX

1: SharpsAway® Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX

1: SharpsAway® II Locking Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX
- 1: Drape: 41" x 55" with 4" fenestration Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX

1: Towel: 13" x 18" (33 cm x 46 cm) Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX

1: Needle Holder Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX

1: Safety Scalpel: #11 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX

1: Medication Label: 1% Lidocaine Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX

4: Gauze Pad: 2" x 2" (5 cm x 5 cm) Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX

4: Gauze Pad: 4" x 4" (10 cm x 10 cm) Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX

1: Dressing: Tegaderm® 3 10 cm x 12 cm Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX

1: Surgical Apparel: Impervious Gown Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX

1: Surgical Apparel: Bouffant Cap Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX

1: Surgical Apparel: Mask with Eye Shield Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX

1: Suture: 3-0 Silk with Curved Needle Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX

²A registered trademark of CareFusion or one of its subsidiaries.
³A registered trademark of 3M Company.

Catheter Priming Volume: 0.24 mL Priming volume is approximate.

Not made with natural rubber latex.
Store between 20 – 25°C (68 – 77°F).
Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

