



SAC-01618-X1A

(IPNXXXXXX)

18^{Ga.}

16^{cm} catheter
length

.025^{inch dia.}
spring-wire
guide

ARROW®

ErgoPack®



Arterial Catheterization

Contents:

- | | | |
|--|---|---|
| 1: Catheter: 18 Ga. x 16 cm Radiopaque PEBA Origin: XX XXXX XX XXX
XXXXXXXXXXXXXXXXXXXX | 1: Maximal Barrier Drape™ with 4" fenestration Origin: XX XXXX XX XXX
XXXXXXXXXXXXXXXXXXXX | XXXXXXXXXXXX |
| 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on
One End - "J" Tip on Other) with Arrow Advancer Origin: XX XXXX XX XXX
XXXXXXXXXXXXXXXXXXXX | 1: Towel: 13" x 18" (33 cm x 46 cm) Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX | 1: HemoHopper® Fluid Receptacle Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX |
| 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW Origin: XX XXXX XX XXX
XXXXXXXXXXXXXXXXXXXX | 1: Needle Holder Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | ¹ A trademark of Becton, Dickinson and Company. |
| 1: Injection Needle: SafetyGlide™ ¹ 25 Ga. x 1" (2.54 cm) Origin: XX XXXX XX
XXX XXXXXXXX XXXXXXXXXXXX | 1: Safety Scalpel: #11 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | ² A registered trademark of CareFusion or one of its subsidiaries. |
| 1: Injection Needle: SafetyGlide™ ¹ 23 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX
XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Checklist Sheet Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | ³ A registered trademark of 3M Company. |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX
XXXXXXXXXXXXXXXXXXXX | 1: Medication Label: 1% Lidocaine Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX | ⁴ A registered trademark of Johnson & Johnson Corporation. |
| 1: Syringe: 3 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX | |
| 1: Syringe: 5 mL Luer-Slip Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX | |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep Origin: XX XXXX XX XXX
XXXXXXXXXXXXXXXXXXXX | 1: Dressing: Tegaderm® ³ 10 cm x 12 cm Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX | |
| 2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® ² One-Step Solution with
Hi-Lite Orange™ Tint Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | 1: Surgical Apparel: Impervious Gown Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX | |
| 1: SharpsAway® Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX | 1: Dressing: BIOPATCH® ⁴ Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX | |
| 1: SharpsAway® II Locking Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX | 1: Surgical Apparel: Bouffant Cap Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX | |
| | 1: Surgical Apparel: Mask with Eye Shield Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX | |
| | 1: Suture: 3-0 Silk with Curved Needle Origin: XX XXXX XX XXX XXXXXXXX | |

Catheter Priming Volume: 0.24 mL Priming volume is approximate.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

