

All Products with Component: Catheter: 20 Ga. x 1-1/2"
(3.81 cm) Radiopaque Polyurethane over 21 Ga. TW
Introducer Needle with Integral .018" (0.46 mm) dia.
Spring-Wire Guide

**QuickFlash® Radial Artery Catheterization**

RA CATH KIT: 20 GA X 1-1/2"

SKU / ARTICLE #: AK-04220**QuickFlash® Radial Artery Catheterization Set**

RA CATH SET: 20 GA X 1 1/2IN

SKU / ARTICLE #: RA-04220-W