

All Products with Component: Three Lumen TaperFree<sup>®</sup>  
Catheter: 6 Fr. (2.14 mm OD) x 40 cm, Pressure  
Injectable, T-Port Connector, Blue FlexTip<sup>®</sup>,  
Contamination Guard and Placement Wire

NO IMAGE  
AVAILABLE

Pressure Injectable Three-Lumen PICC

PI PICC KIT: 3-L 6 FR X 40 CM

**SKU / ARTICLE #: CDC-34063-HPK1A**