



## Trauma

### TRAUMA KIT: 8.5 FR

SKU / Article #: AK-05801

### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 10**



*Representative photo(s): actual product specifications may vary from that shown in the photo*

Trauma SKU / Article #: AK-05801	
Quantity	Component Description
1	Sheath: 8.5 Fr. x 3-1/2" (8.89 cm) Arrow-Flex® Radiopaque FEP with Tissue Dilator
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Pressure Transduction Probe
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Syringe: 5 mL Luer-Slip
1	Syringe: 30 mL Luer-Lock
1	3 mL Applicator 2% CHG and 70% IPA ChloroPrep® <sup>1</sup> One-Step Solution with Hi-Lite Orange™ Tint
1	SharpsAway® II Locking Disposal Cup
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
1	Scalpel: #11
1	Dressing: OpSite® <sup>2</sup> 14 cm x 10 cm
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup> A registered trademark of CareFusion or one of its subsidiaries.

<sup>2</sup> A registered trademark of T.J. Smith & Nephew, Limited.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.