



Peritoneal Lavage

PERITONEAL LAVAGE CATH KIT:8FR

SKU / Article #: AK-09001

Considerations




Not made with natural rubber latex.
 ±Precaution: Contains Phthalate: DEHP

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

Peritoneal Lavage SKU / Article #: AK-09001	
Quantity	Component Description
1	Peritoneal Lavage Catheter: 8 Fr. x 9-1/4" (24 cm) Radiopaque FEP with Side Holes
1	Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
2	Syringe: 5 mL Luer-Slip
1	Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
1	Scalpel: #11
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Tubing Set±: Fluid Administration

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.