



Arrowgard Blue[®] PSI for use with 7.5 - 8 Fr. Catheters

PSI KIT: 9 FR

SKU / Article #: AK-29907-CDC

Highlights

Antimicrobial: Arrowgard Blue*
Sharps Safety

Considerations

*Chlorhexidine
Not made with natural rubber latex.
Includes Lidocaine
±Precaution: Contains Phthalate: DEHP

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowgard Blue [®] PSI for use with 7.5 - 8 Fr. Catheters SKU / Article #: AK-29907-CDC	
Quantity	Component Description
1	Single-Lumen Infusion Catheter, SLIC [®] : 7 Fr. x 16 cm
1	Arrowgard Blue [®] Arrow-Flex [®] Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Pressure Transduction Probe
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow [®] Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe





Arrowgard Blue PSI for use with 7.5 - 8 Fr. Catheters

SKU / Article #: AK-29907-CDC

Quantity	Component Description
1	Syringe: 3 mL Luer-Lock
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA ChloraPrep ^{®1} One-Step Solution with Hi-Lite Orange [™] Tint
1	Obturator: 8 Fr.
1	Cath-Gard [®] Catheter Contamination Shield: 80 cm with TwistLock [™] Adapter
1	Dust Cap: Non-Vented
1	SharpsAway [®] II Locking Disposal Cup
1	SharpsAway [®] Disposal Cup
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
1	Drape: 24" x 36" (61 cm x 91 cm)
1	Drape: 41" x 55" with 4" fenestration
1	Needle Holder
1	Safety Scalpel: #11
1	Stopcock: 4-way High-Flow
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Silk with Curved Needle
1	Suture: 3-0 Braided Silk with Straight Needle
1	Tubing±: 20" (50 cm) Extension

¹ A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.