



Arterial Catheterization Kit

ARTERIAL CATH KIT: 20GA X 5CM

SKU / Article #: ASK-00520-WHO1

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Arterial Catheterization Kit	
SKU / Article #: ASK-00520-WHO1	
Quantity	Component Description
1	Catheter: 20 Ga. x 5 cm Radiopaque PEBA
1	Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Injection Needle: Eclipse™ ¹ 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Syringe: 5 mL Luer-Slip
1	Dust Cap: Non-Vented
1	SharpsAway® II Locking Disposal Cup
1	Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
3	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® ² CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)

¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of 3M Company.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.