

PSI Kit for use with 7.5 - 8 Fr. Catheters

PSI KIT: 8.5 FR

SKU / Article #: ASK-09803-FM1

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

NO IMAGE
AVAILABLE

Representative photo(s): actual product specifications may vary from that shown in the photo

PSI Kit for use with 7.5 - 8 Fr. Catheters


SKU / Article #: ASK-09803-FM1

Quantity	Component Description
1	Single-Lumen Infusion Catheter, SLIC®: 7 Fr. x 16 cm
1	Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Injection Needle: Eclipse™ ¹ 23 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: Eclipse™ ¹ 25 Ga. x 1" (2.54 cm)
1	Syringe: 3 mL Luer-Lock
2	Syringe: 5 mL Luer-Slip
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
1	SharpsAway® Disposal Cup
1	Maximal Barrier Drape™ with 4" fenestration
1	Safety Scalpel: #11
1	Dressing: Tegaderm® ² 4" x 4-3/4" (10 cm x 12 cm)
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Transducer Cover: 14 cm x 147 cm
1	Suture: 3-0 Braided Silk with Straight Needle

¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of 3M Company.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

© 2024 Teleflex Incorporated. All Rights Reserved.