



Arrowgard Blue® Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions

HEMODIALYSIS KIT: 2L 12 FR X 16 CM

SKU / Article #: CA-22122-F

Highlights

Antimicrobial: Arrowgard Blue*

Considerations

*Chlorhexidine

Not made with natural rubber latex.

Sales Quantity/Case: 5




Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowgard Blue® Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions

SKU / Article #: CA-22122-F

Quantity	Component Description
1	Two-Lumen Catheter: 12 Fr. (4.06 mm OD) x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowgard® Antimicrobial Surface Treatment, Extension Line Clamps
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Pressure Transduction Probe
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
2	Dust Cap: Non-Vented
1	SharpsAway® Disposal Cup
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
1	Scalpel: #11
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.