

*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Arrowg+ard Blue® Large-Bore Three-Lumen CVC Kit for High Volume Infusions

**CVC KIT: 3 LUMEN 12 FR X 20 CM**

SKU / Article #: CA-25123-F

### Highlights

Antimicrobial: Arrowg+ard Blue\*

### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.



**Sales Quantity/Case: 5**

# Arrow+ard Blue® Large-Bore Three-Lumen CVC Kit for High Volume Infusions

SKU / Article #: CA-25123-F

Quantity	Component Description
1	Arrow+ard Blue® Three-Lumen Catheter: 12 Fr. (4.1 mm OD) x 20 cm
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Pressure Transduction Probe
1	Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
3	Dust Cap: Non-Vented
1	SharpsAway® Disposal Cup
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
1	Scalpel: #11
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for

Use available under the Documents section on this page, where an IFU is available.

© 2025 Teleflex Incorporated. All Rights Reserved.