

NO IMAGE  
AVAILABLE

## Arterial Catheterization Kit

### ARTERIAL CATH KIT: 20 GA

SKU / Article #: CAN-00820-NSH

#### Highlights

---

Sharps Safety

#### Considerations

---

Not made with natural rubber latex.

±Precaution: Contains Phthalate: DEHP

**Sales Quantity/Case: 5**

---

*Representative photo(s): actual product specifications may vary from that shown in the photo*

# Arterial Catheterization Kit






SKU / Article #: CAN-00820-NSH

Quantity	Component Description
1	Catheter: 20 Ga. x 8 cm Radiopaque Polyurethane
1	Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Injection Needle: SafetyGlide™ <sup>1</sup> 18 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)
1	Syringe: 10 mL Luer-Lock
1	Dust Cap: Non-Vented
1	SharpsAway® Disposal Cup
1	SharpsAway® II Locking Disposal Cup
1	Drape: 41" x 55" with 4" fenestration
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Transducer Cover: 14 cm x 147 cm
1	Dressing: STATLOCK® <sup>2</sup> PICC Plus Catheter Stabilization Device
1	Tubing±: 24" (61 cm) Pressure Monitoring

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of C. R. Bard, Inc.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

