

NO IMAGE  
AVAILABLE

## Arterial Catheterization Kit

ARTERIAL CATH KIT: 20 GA X 12 CM

SKU / Article #: CAN-01220-IUC

### Highlights

---

Sharps Safety

### Considerations

---

Not made with natural rubber latex.

±Precaution: Contains Phthalate: DEHP

**Sales Quantity/Case: 5**

---

*Representative photo(s): actual product specifications may vary from that shown in the photo*




Arterial Catheterization Kit  
SKU / Article #: CAN-01220-IUC

Quantity	Component Description
1	Catheter: 20 Ga. x 12 cm Radiopaque PEBA
1	Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Injection Needle: SafetyGlide™ <sup>1</sup> 18 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)
1	Syringe: 10 mL Luer-Lock
1	Syringe: 5 mL Luer-Slip
1	Dust Cap: Non-Vented
1	SharpsAway® Disposal Cup
1	SharpsAway® II Locking Disposal Cup
1	Drape: 41" x 55" with 4" fenestration
1	Needle Holder
1	Safety Scalpel: #11
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® <sup>2</sup> 8.5 cm x 10.5 cm
1	Transducer Cover: 14 cm x 147 cm
1	Suture: 3-0 Silk with Curved Needle
1	Tubing±: 24" (61 cm) Pressure Monitoring

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of 3M Company.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.