

## Pressure Injectable Two-Lumen PICC

### PI PICC KIT: 2-L 5 FR X 50 CM

SKU / Article #: CDA-35052-HPK1A

#### Highlights

Sharps Safety  
CE Mark

#### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 5**



*Representative photo(s): actual product specifications may vary from that shown in the photo*

### Pressure Injectable Two-Lumen PICC

SKU / Article #: CDA-35052-HPK1A

Quantity	Component Description
1	Two Lumen TaperFree® Catheter: 5 Fr. (1.81 mm OD) x 50 cm, Pressure Injectable, T-Port Connector, Blue FlexTip®, Contamination Guard and Placement Wire
1	GlideThru™ Peel-Away Sheath: 5 Fr. x 2-3/4" (7 cm) Radiopaque over 5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: SafetyGlide™1 25 Ga. x 5/8" (1.60 cm) RW
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup

## Pressure Injectable Two-Lumen PICC






SKU / Article #: CDA-35052-HPK1A

Quantity	Component Description
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Patient ID Card
1	Sterile Procedure Sign
1	Chart Sticker
1	Patient Information Booklet
1	Checklist/CLIP Sheet
2	Paper Tape Measure
1	Tourniquet
5	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Surgical Apparel: Mask with Eye Shield
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask
1	Tape: Steri-Strip® <sup>2</sup>
1	HemoHopper® Fluid Receptacle

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of 3M Company.

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.