



Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC pre-loaded with Arrow® VPS TipTracker™ Stylet

PI PICC KIT: 1L 4.5FR 55CM AGBA TIPTRACK

SKU / Article #: CDA-45541-TTS

Highlights

- Antimicrobial: Arrowg+ard Blue Advance*
- Antithrombogenic: Arrowg+ard Blue Advance*
- Sharps Safety
- CE Mark
- Pre-loaded with VPS Stylet

Considerations

- *Chlorhexidine
- Not made with natural rubber latex.

Sales Quantity/Case: 3

Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC pre-loaded with Arrow® VPS TipTracker™ Stylet

SKU / Article #: CDA-45541-TTS

Quantity	Component Description
1	One Lumen TaperFree® Catheter with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 4.5 Fr. (1.58 mm OD) x 55 cm, Pressure Injectable, Blue FlexTip®, Arrow® TipTracker™ Stylet and T-Port Connector
1	GlideThru™ Peel-Away Sheath: 4.5 Fr. x 2-3/4" (7 cm) Radiopaque over 4.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: Eclipse™ ² 25 Ga. x 1" (2.54 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock

Pressure Injectable Arrowg+ard Blue Advance One-Lumen PICC pre-loaded with Arrow VPS
TipTracker Stylet

SKU / Article #: CDA-45541-TTS

Quantity	Component Description
1	Syringe: 3 mL Luer-Lock
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Patient ID Card
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
1	Foam Electrodes (3 per pack)
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm®3 15.5 cm x 10 cm
1	Surgical Apparel: Impervious Gown
1	Dressing: STATLOCK®4 Catheter Stabilization Device
1	Surgical Apparel: Mask with Eye Shield
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask
1	Equipment Cover: 13 cm x 244 cm (5" x 96")
1	Tape: Steri-Strip®3
1	HemoHopper® Fluid Receptacle

Pressure Injectable Arrowgard Blue Advance One-Lumen PICC pre-loaded with Arrow VPS TipTracker Stylet
 SKU / Article #: CDA-45541-TTS

Quantity	Component Description
1	T-Piece Cover

¹ Licensed under US Patent No. 7,329,412.

² A trademark of Becton, Dickinson and Company.

³ A registered trademark of 3M Company.

⁴ A registered trademark of C. R. Bard, Inc.

Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.