

Pressure Injectable Two-Lumen PICC pre-loaded with Arrow® VPS® Stylet PI PICC:2L 5FR X 50CM PRELOAD VPS STYLET

SKU / Article #: CDC-35052-VPS

Highlights

Sharps Safety
Pre-loaded with VPS Stylet

Considerations

Not made with natural rubber latex.
Includes Lidocaine

Sales Quantity/Case: 3

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Two-Lumen PICC pre-loaded with Arrow® VPS® Stylet SKU / Article #: CDC-35052-VPS	
Quantity	Component Description
1	Two Lumen TaperFree® Catheter: 5 Fr. (1.81 mm OD) x 50 cm, Pressure Injectable, Blue FlexTip®, Contamination Guard, Arrow® VPS® Stylet and Tuohy-Borst Adapter
1	GlideThru™ Peel-Away Sheath: 5 Fr. x 2-3/4" (7 cm) Radiopaque over 5 Fr. Dilator
1	Placement Wire Assembly
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: SafetyGlide™1 25 Ga. x 5/8" (1.60 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock

Pressure Injectable Two-Lumen PICC pre-loaded with Arrow VPS Stylet

SKU / Article #: CDC-35052-VPS


Quantity	Component Description
4	10 mL Luer-Lock Pre-Filled Saline Syringe
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep ^{®2} One-Step Solution
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep ^{®2} One-Step Solution with Hi-Lite Orange [™] Tint
2	Dust Cap: Non-Vented
1	SecondSite [™] Adjustable Hub: Fastener
1	SecondSite [™] Adjustable Hub: Catheter Clamp
1	SharpsAway [®] II Locking Disposal Cup
1	SharpsAway [®] Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape [™] with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Patient ID Card
1	Chart Sticker
1	Patient Information Booklet
1	Medication Label: 1% Lidocaine
1	Checklist/CLIP Sheet
2	Paper Tape Measure
1	Tourniquet
1	Foam Electrodes (3 per pack)
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm

Pressure Injectable Two-Lumen PICC pre-loaded with Arrow VPS Stylet
 SKU / Article #: CDC-35052-VPS

Quantity	Component Description
1	Dressing: BIOPATCH ^{®3}
1	Surgical Apparel: Mask with Eye Shield
1	Dressing: STATLOCK ^{®4} PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask
1	Tape: Steri-Strip ^{®5}
1	HemoHopper [®] Fluid Receptacle

- ¹ A trademark of Becton, Dickinson and Company.
- ² A registered trademark of CareFusion or one of its subsidiaries.
- ³ A registered trademark of Johnson & Johnson Corporation.
- ⁴ A registered trademark of C. R. Bard, Inc.
- ⁵ A registered trademark of 3M Company.

Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.