

Three-Lumen CVC

CVC SET: 3-LUMEN 8.5FR X 11CM

SKU / Article #: CS-10853

Highlights

CE Mark

Considerations




Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Three-Lumen CVC	
SKU / Article #: CS-10853	
Quantity	Component Description
1	Three-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 11 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
3	Dust Cap: Non-Vented
3	Clamp: Extension Line
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.