

Pediatric Arrowgard Blue® Two-Lumen CVC

CVC SET: 2-LUMEN 5 FR X 8 CM

SKU / Article #: CS-25502

Highlights

Antimicrobial: Arrowgard Blue*
CE Mark

Considerations




*Chlorhexidine
Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Pediatric Arrowgard Blue® Two-Lumen CVC	
SKU / Article #: CS-25502	
Quantity	Component Description
1	Arrowgard Blue® Two-Lumen Catheter: 5 Fr. (1.8 mm OD) x 8 cm
1	Spring-Wire Guide: .021" (0.53 mm) dia. x 17-3/4" (45 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other)
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
1	Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.