



Three-Lumen CVC

CVC SET: 3-LUMEN 8.5FR X 20CM

SKU / Article #: DE-15853-S

Highlights

Sharps Safety
CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Three-Lumen CVC	
SKU / Article #: DE-15853-S	
Quantity	Component Description
1	Three-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 20 cm
1	Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
3	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

