

## Pressure Injectable Arrowg+ard Blue Advance® Three-Lumen PICC pre-loaded with Arrow® VPS Precision™ Stylet

AGBA PICC G4 Stylet: 3L 6FR x 55CM TCG

SKU / Article #: DLX-45563-VPSC

### Highlights

- Antimicrobial: Arrowg+ard Blue Advance\*
- Antithrombogenic: Arrowg+ard Blue Advance\*
- Sharps Safety
- Pre-loaded with VPS Stylet

### Considerations

- \*Chlorhexidine
- Not made with natural rubber latex.
- Includes Lidocaine
- ###MR Safe

**Sales Quantity/Case: 3**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Pressure Injectable Arrowg+ard Blue Advance® Three-Lumen PICC pre-loaded with Arrow® VPS Precision™ Stylet

SKU / Article #: DLX-45563-VPSC

Quantity	Component Description
1	Three Lumen TaperFree® Catheter### with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection <sup>1</sup> : 6 Fr. (2.11 mm OD) x 55 cm, Pressure Injectable, Contamination Guard, Arrow® VPS Precision™ Stylet, Tuohy-Borst Adapter, Stiffening Stylet (0.36 mm dia x 75 cm with straight distal tip) and T-Port Connector
1	GlideThru™ Peel-Away Sheath: 6 Fr. x 2-3/4" (7 cm) Radiopaque over 6 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: SafetyGlide™ <sup>2</sup> 25 Ga. x 5/8" (1.60 cm) RW
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW

Pressure Injectable Arrow+ard Blue Advance Three-Lumen PICC pre-loaded with Arrow VPS Precision Stylet  
 SKU / Article #: DLX-45563-VPSC

Quantity	Component Description
1	Syringe: 5 mL Luer-Lock
4	10 mL Luer-Lock Pre-Filled Saline Syringe
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep <sup>®3</sup> One-Step Solution
2	3 mL Applicator 2% CHG and 70% IPA Chloraprep <sup>®3</sup> One-Step Solution with Hi-Lite Orange <sup>™</sup> Tint
3	Dust Cap: Non-Vented
1	SharpsAway <sup>®</sup> II Locking Disposal Cup
1	SharpsAway <sup>®</sup> Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape <sup>™</sup> with two 13 cm (5") x 15 cm (6") fenestrations
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Patient ID Card
1	Patient Information Booklet
1	Medication Label: 1% Lidocaine
2	Paper Tape Measure
1	Tourniquet
1	Foam Electrodes (3 per pack)
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: Tegaderm <sup>®4</sup> CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)
1	Transducer Cover with Gel and Rubber Bands
1	Surgical Apparel: Mask with Eye Shield
1	Dressing: STATLOCK <sup>®5</sup> PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap

Pressure Injectable Arrowgard Blue Advance Three-Lumen PICC pre-loaded with Arrow VPS Precision Stylet  
 SKU / Article #: DLX-45563-VPSC

Quantity	Component Description
1	Surgical Apparel: Mask
1	Tape: Steri-Strip® <sup>4</sup>

<sup>1</sup> Licensed under US Patent No. 7,329,412 and 11,219,706.


<sup>2</sup> A trademark of Becton, Dickinson and Company.

<sup>3</sup> A registered trademark of CareFusion or one of its subsidiaries.

<sup>4</sup> A registered trademark of 3M Company.

<sup>5</sup> A registered trademark of C. R. Bard, Inc.

## Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.