

Home / Midline and Peripheral IV / **EDC-00820-HB**



Arrow™ Endurance™ Extended Dwell Peripheral Catheter System

EXT DWELL CATH DEVICE 20 GA X 8 CM

SKU / Article #: EDC-00820-HB

Highlights


Sharps Safety

Considerations

Not made with natural rubber latex.

††MR Safe

Sales Quantity/Case: 5




Representative photo(s): actual product specifications may vary from that shown in the photo

Arrow™ Endurance™ Extended Dwell Peripheral Catheter System
SKU / Article #: EDC-00820-HB

Quantity	Component Description
1	Catheter System: 20 Ga. (1.10 mm) x 80 mm (3.15") Arrow™ Endurance™ over 22 Ga. (0.8 mm) RW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
1	10 mL Luer-Lock Pre-Filled Saline Syringe
1	3 mL Applicator 2% CHG and 70% IPA ChloraPrep® ¹ One-Step Solution with Hi-Lite Orange™ Tint
1	Arm Drape: 36" x 46" with 3" x 5" Fenestration
1	Label: Venous/Arterial Placement
1	Paper Tape Measure
1	Tourniquet
2	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Dressing: BIOPATCH® ²
1	Transducer Cover with Gel and Rubber Bands
1	Dressing: STATLOCK® ³ PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask

¹ A registered trademark of CareFusion or one of its subsidiaries.
² A registered trademark of Johnson & Johnson Corporation.
³ A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

