

## Peripheral Emergency Infusion Device - EID™

### EMERGENCY INFUSION DEVICE SET: 6 FR

SKU / Article #: EI-04060

### Highlights

CE Mark

### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 10**


*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Peripheral Emergency Infusion Device - EID™

SKU / Article #: EI-04060

Quantity	Component Description
1	Sheath: 6 Fr. x 1-1/2" (3.81 cm) Radiopaque with Dilator over 20 Ga. TW Needle
1	Spring-Wire Guide/Tube Assembly: .025" (0.64 mm) dia. x 6-7/8" (17.5 cm)
1	Scalpel: #11

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.