



Emergency Infusion Device - EID™

EMERGENCY INFUSION DEVICE SET

SKU / Article #: EI-04080

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10


Representative photo(s): actual product specifications may vary from that shown in the photo

Emergency Infusion Device - EID™

SKU / Article #: EI-04080

Quantity	Component Description
1	EID™ Emergency Infusion Device: 8.5 Fr. x 3-1/2" (8.9 cm) ArrowFlex® Radiopaque FEP Sheath/Dilator over 17 Ga. x 5-1/2" (14 cm) TW Needle, .035" (0.89 mm) dia. x 13-3/16" (33.5 cm) Spring-Wire Guide, 5 mL Luer-Lock Syringe and Sterility Tubing
1	Scalpel: #11

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

© 2024 Teleflex Incorporated. All Rights Reserved.