

One-Lumen CVC

CVC SET: 22 GA X 4" (10 CM)

SKU / Article #: ES-04522


Highlights

CE Mark

Considerations

Not made with natural rubber latex.




Sales Quantity/Case: 10



Representative photo(s): actual product specifications may vary from that shown in the photo

One-Lumen CVC	
SKU / Article #: ES-04522	
Quantity	Component Description
1	One-Lumen Catheter: 22 Ga. (0.9 mm OD) x 10 cm
1	Spring-Wire Guide: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 4 Fr. (1.35 mm) x 7.6 cm
1	Tissue Dilator: 4 Fr. (1.35 mm) x 4.1 cm
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.