



PSI Kit for use with 7 - 7.5 Fr. Catheters

PSI KIT: 8 FR

SKU / Article #: ES-09807

Considerations

Caution: This product contains dry natural rubber.
Includes Lidocaine

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo



PSI Kit for use with 7 - 7.5 Fr. Catheters

SKU / Article #: ES-09807

Quantity	Component Description
1	Arrow-Flex® Sheath: 8 Fr. (2.5 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Syringe: 3 mL Luer-Lock
1	Syringe: 5 mL Luer-Slip
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep® ¹ One-Step Solution with Hi-Lite Orange™ Tint
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with Tuohy-Borst Adapter, Sterile Tape
1	SharpsAway® Disposal Cup
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration, with adhesive
2	Towel
1	Scalpel: #11
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 1-0 Silk with Straight Needle

¹ A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.