

Arrow® Midline Kit

MIDLINE CATHETERIZATION KIT

SKU / Article #: EU-02031-ML

Highlights

Sharps Safety
CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5



Representative photo(s): actual product specifications may vary from that shown in the photo

Arrow® Midline Kit	
SKU / Article #: EU-02031-ML	
Quantity	Component Description
1	Midline Catheter, Trimmable: 3 Fr. (1.0 mm OD) x 20 cm, T-Port Connector, Contamination Guard and Placement Wire
1	GlideThru™ Peel-Away Sheath: 3 Fr. x 2-3/4" (7 cm) Radiopaque over 3 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	SharpsAway® II Locking Disposal Cup
1	Catheter Trimmer
1	Safety Scalpel: #11
1	Paper Tape Measure
1	Tourniquet
1	Dressing: GRIP-LOK® ¹ Catheter Fixation Device
1	Valve: MicroCLAVE™ ² Neutral Displacement Connector

¹ A registered trademark of TIDI.

² A trademark of ICU Medical, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.