

## Two-Lumen CVC

### CVC SET: 2-LUMEN 4FR X 13CM

SKU / Article #: EU-14402-N

### Highlights

Sharps Safety  
CE Mark

### Considerations


Not made with natural rubber latex.

**Sales Quantity/Case: 5**

Two-Lumen CVC	
SKU / Article #: EU-14402-N	
Quantity	Component Description
1	Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 13 cm
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
1	Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

*Representative photo(s): actual product specifications may vary from that shown in the photo*

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.