

*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Pressure Injectable Arrowg+ard Blue Plus® Four-Lumen CVC

CVC SET: 4-LUMEN 8.5FR X 16CM

SKU / Article #: EU-42854-HPS

### Highlights

- Antimicrobial: Arrowg+ard Blue Plus\*
- Sharps Safety
- CE Mark

### Considerations

- \*Chlorhexidine
- Not made with natural rubber latex.

**Sales Quantity/Case: 10**




### Pressure Injectable Arrowg+ard Blue Plus® Four-Lumen CVC

SKU / Article #: EU-42854-HPS

Quantity	Component Description
1	Arrowg+ard Blue Plus® <sup>1</sup> Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 16 cm, Pressure Injectable
1	Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 9 Fr. x 10.2 cm
4	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

<sup>1</sup> Licensed under US Patent No. 7,329,412.

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.