

Arterial Catheterization Set

ARTERIAL CATH SET: 20 GA X 3-1/16 IN

SKU / Article #: GH-04125

Highlights

CE Mark

Considerations


Not made with natural rubber latex.

Sales Quantity/Case: 50

Representative photo(s): actual product specifications may vary from that shown in the photo

Arterial Catheterization Set SKU / Article #: GH-04125	
Quantity	Component Description
1	Catheter: 20 Ga. x 3-1/16" (7.78 cm) with Integral Suture Wing
1	Spring-Wire Guide: .025" (0.64 mm) dia. x 9-3/4" (25 cm) (Straight Soft Tip Both Ends)
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

© 2025 Teleflex Incorporated. All Rights Reserved.