

Pressure Injectable One-Lumen PICC

PI PICC: 1L 4FR x 55CM

SKU / Article #: IDLX-35541-HPK

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.
 ###MR Safe

Sales Quantity/Case: 3



Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable One-Lumen PICC SKU / Article #: IDLX-35541-HPK

Quantity	Component Description
1	One Lumen TaperFree® Catheter###: 4 Fr. (1.45 mm OD) x 55 cm, Pressure Injectable, Stiffening Stylet (0.36 mm dia. x 75 cm with straight distal tip) and T-Port Connector
1	GlideThru™ Peel-Away Sheath: 4 Fr. x 2-3/4" (7 cm) Radiopaque over 4 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Injection Needle: Eclipse™1 25 Ga. x 1" (2.54 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp

Pressure Injectable One-Lumen PICC

SKU / Article #: IDLX-35541-HPK

Quantity	Component Description
1	SecondSite™ Adjustable Hub: Fastener
1	SharpsAway® Disposal Cup
1	SharpsAway® II Locking Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Patient ID Card
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
5	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Transducer Cover with Gel and Rubber Bands
1	Surgical Apparel: Mask with Eye Shield
1	Dressing: STATLOCK® ² PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask
1	Tape: Steri-Strip® ³
1	Valve: MicroCLAVE™ ⁴ Neutral Displacement Connector







¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of C. R. Bard, Inc.

³ A registered trademark of 3M Company.

⁴ A trademark of ICU Medical, Inc.

Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.