



ARROW[®] GlideThru[™] Peel-Away Sheath over Dilator Dilator accepts .018" dia. or smaller Spring-Wire Guide. Sheath accepts 6 Fr. or smaller catheter.

6FR X 10CM GLIDETHRU SHEATH/DILATOR

SKU / Article #: PL-01060X

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

ARROW[®] GlideThru[™] Peel-Away Sheath over Dilator Dilator accepts .018" dia. or smaller Spring-Wire Guide. Sheath accepts 6 Fr. or smaller catheter.

SKU / Article #: PL-01060X

Quantity	Component Description
1	GlideThru [™] Peel-Away Sheath: 6 Fr. x 4" (10 cm) Radiopaque over 6 Fr. Dilator

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

