

Pressure Injectable Two-Lumen Midline

PI MIDLINE 2-L: 5FR X 20CM

SKU / Article #: PR-32052-BAS

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.

††MR Safe

Sales Quantity/Case: 5



Representative photo(s): actual product specifications may vary from that shown in the photo


Pressure Injectable Two-Lumen Midline

SKU / Article #: PR-32052-BAS

| Quantity | Component Description |
|----------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | Two-Lumen TaperFree® Catheter†††: 5 Fr.(1.81 mm OD) x 20 cm Pressure Injectable Radiopaque Polyurethane with Extension Line Clamps and Contamination Guard |
| 1 | GlideThru™ Peel-Away Sheath: 5 Fr. x 2-3/4" (7 cm) Radiopaque over 5 Fr. Dilator |
| 1 | Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) |
| 1 | Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW |
| 1 | Syringe: 10 mL Luer-Lock |
| 2 | Dust Cap: Non-Vented |
| 1 | SharpsAway® II Locking Disposal Cup |
| 1 | SharpsAway® Disposal Cup |
| 1 | Catheter Trimmer |
| 1 | Safety Scalpel: #11 |
| 1 | Overbed Sign |
| 1 | Patient Information Booklet |
| 1 | Paper Tape Measure |
| 1 | Dressing: STATLOCK® ¹ PICC Plus Catheter Stabilization Device |

¹ A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

