

4Fr. x 50cm Pressure Injectable Multi-Lumen ARROW® PICC Set with 130cm Stainless Steel Marked Wire

PI PICC 1-LUMEN: 4FR X 50CM W/ 130CM SS SWG

SKU / Article #: PR-35041-HPSL

## Contraindications & Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

*Representative photo(s): actual product specifications may vary from that shown in the photo*



## 4Fr. x 50cm Pressure Injectable Multi-Lumen ARROW® PICC Set with 130cm Stainless Steel Marked Wire

SKU / Article #: PR-35041-HPSL

Quantity	Component Description
1	Indwelling TaperFree™ Catheter: 4 Fr. x 50 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip® and Extension Line Clamp
1	GlideThru™ Peel-Away Sheath: 4 Fr. x 4" (10 cm) Radiopaque over 4 Fr. Dilator
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)
1	Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 51-3/16" (130 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Lock
1	Skin Protectant Prep Pad (non-sterile solution)
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	Flow Rate Information Card with Injection Log
1	Patient ID Card
1	Chart Sticker
1	Patient Information Booklet
1	Checklist/CLIP Sheet
1	Paper Tape Measure
1	Dressing: STATLOCK® <sup>1</sup> Catheter Stabilization Device

<sup>1</sup> A registered trademark of C. R. Bard, Inc.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.