

*Representative photo(s): actual product specifications may vary from that shown in the photo*

## MAC™ Two-Lumen Central Venous Access Set with Integral Hemostasis Valve for use with 7 - 7.5 Fr. Catheters

### MULTI-LUMEN/PSI SET

SKU / Article #: SI-11242

### Highlights

---

CE Mark

### Considerations

---

Not made with natural rubber latex.

**Sales Quantity/Case: 5**



---

# MAC™ Two-Lumen Central Venous Access Set with Integral Hemostasis Valve for use with 7 - 7.5 Fr. Catheters

SKU / Article #: SI-11242

Quantity	Component Description
1	Two-Lumen Central Venous Access Device: 4-1/2" (11.5 cm) Radiopaque Polyurethane with 9 Fr. Distal Lumen, Integral Hemostasis Valve, Extension Line Clamps
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with Tuohy-Borst Adapter, Sterile Tape
1	Dust Cap: Vented
3	Dust Cap: Non-Vented
1	Stopcock: 4-way High-Flow
3	Gauze Pad: 4" x 4" (10 cm x 10 cm)

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

