

PSI Kit for use with 7 - 7.5 Fr. Catheters

PSI KIT: 8.5 FR

SKU / Article #: SP-09801

Highlights

CE Mark

Considerations

Caution: This product contains dry natural rubber.

Sales Quantity/Case: 10

NO IMAGE
AVAILABLE



Representative photo(s): actual product specifications may vary from that shown in the photo

PSI Kit for use with 7 - 7.5 Fr. Catheters

SKU / Article #: SP-09801

Quantity	Component Description
1	Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 3 mL Luer-Lock
2	Syringe: 5 mL Luer-Slip
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with Tuohy-Borst Adapter, Sterile Tape
1	Dust Cap: Non-Vented
1	SharpsAway® Disposal Cup
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
1	Scalpel: #11
1	Stopcock: 3-way High-Flow
2	Prep Sponge Swab
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle
1	Tray: Prep

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

© 2024 Teleflex Incorporated. All Rights Reserved.