



Representative photo(s): actual product specifications may vary from that shown in the photo

# Percutaneous Sheath Introducer for use with 7 - 8 Fr. Catheters

PSI COMPONENT: 9 FR

SKU / Article #: SS-09903-S

## Highlights

CE Mark


## Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 25

Percutaneous Sheath Introducer for use with 7 - 8 Fr. Catheters SKU / Article #: SS-09903-S	
Quantity	Component Description
1	Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 7 cm and Tissue Dilator
1	Obturator: 8 Fr.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

