

NO IMAGE
AVAILABLE

Femoral Arterial Line Catheterization

ACCESS TRAY

SKU / Article #: UM-04018

Considerations


Not made with natural rubber latex.
Includes Lidocaine

Sales Quantity/Case: 10

*Representative photo(s): actual product
specifications may vary from that shown in the
photo*

Femoral Arterial Line Catheterization SKU / Article #: UM-04018	
Quantity	Component Description
1	Indwelling Catheter: 18 Ga. x 16 cm with Integral Suture Wing
1	Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
1	Introducer Needle: 18 Ga. x 3" (7.62 cm) RW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
2	Syringe: 5 mL Luer-Slip
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
1	Scalpel: #11
1	Prep Razor
2	Prep Sponge Swab
4	Gauze Pad: 2" x 2" (5 cm x 5 cm)
8	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle
1	Tray: Prep

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for

Use available under the Documents section on this page, where an IFU is available.