



VPS TipTracker™ Stylet Accessory

REF RHY-177-TTSAP

ARROW®

Contents:

- 1: VPS TipTracker™ Stylet with T-Port Connector
0.015" (0.38 mm) dia. x 29" (74 cm)
- 1: Foam Electrodes (3 per pack)
- 1: Equipment Cover: 13 cm x 244 cm (5" x 96")
- 1: T-Piece Cover

Rx only. **Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Sterile. Fluid path components are non-pyrogenic.

- nl** VPS TipTracker stiletaccessoire
- pl** Dodatkowy mandryn TipTracker systemu VPS
- fr** Accessoire stylet VPS TipTracker
- pt** Acessório de estilete VPS TipTracker
- de** VPS TipTracker Mandrin-Zubehörteil
- es** Accesorio de estilete VPS TipTracker
- it** Accessorio con mandrino VPS TipTracker
- sv** VPS TipTracker-mandrängtillbehör



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com