



ARROW ENDURANCE™

Extended Dwell Peripheral Catheter System

22 Ga. Catheter

6 cm catheter length

.010 inch dia. wire guide

Contents:

- 1: Catheter System: 22 Ga. (0.91 mm) x 60 mm (2.36") ARROW ENDURANCE™ over 24 Ga. (0.60 mm) TW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
- 1: Skin Protectant Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint
- 1: Towel: 24" x 36" (61 cm x 91 cm)
- 1: Arm Drape: 36" x 46" with 3" x 5" Fenestration
- 1: Label: Venous/Arterial Placement
- 1: Paper Tape Measure
- 1: Tourniquet
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®² 10 cm x 12 cm
- 1: Dressing: STATLOCK®³ Catheter Stabilization Device
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask

¹A registered trademark of CareFusion or one of its subsidiaries.

²A registered trademark of 3M Company.

³A registered trademark of C. R. Bard, Inc.

Rx only. **Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Priming Volume	Gravity Flow Rate	Pump Flow Rate	MAX Injector Pressure Limit	MAX Pressure Injection Flow Rate per Injectate Viscosity		
				4.7 cP	11.8 cP	26.6 cP
0.37 mL	1422 mL/hr	5109 mL/hr	325 psi	6.0 mL/sec	4.0 mL/sec	2.0 mL/sec

Consult instructions for use for details on the testing performed to generate these values.

Not made with natural rubber latex

Store between 15 - 30°C (59 - 86°F). Avoid freezing and excessive heat above 40°C (104°F).

Arrow International LLC

Subsidiary of Teleflex Incorporated

3015 Carrington Mill Blvd.

Morrisville, NC 27560 USA

Teleflex®

www.teleflex.com

EDC-00622-201B (2021-04)



ARROW ENDURANCE™

Extended Dwell Peripheral Catheter System

22 Ga. Catheter

6 cm catheter length

.010 inch dia. wire guide