



ARROW® PICC Set

REF PR-34063-HPHNM

Contents:

- | | |
|--|---|
| <ul style="list-style-type: none"> 1: Multi-Lumen Indwelling TaperFree™ Catheter: 6 Fr. x 40 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps, and Contamination Guard 1: GlideThru™ Peel-Away Sheath: 6 Fr. x 4" (10 cm) Radiopaque over 6 Fr. Dilator 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends) 1: Guidewire, Marked: .018" (0.46 mm) dia. x 80 cm Hydrophilic-Coated Nitinol with Soft Tip Tungsten Coil (CE 0120) 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW 1: Syringe: 10 mL Luer-Lock 1: Skin Protectant Prep Pad 3: Dust Caps: Non-Vented | <ul style="list-style-type: none"> 1: SecondSite™ Adjustable Hub: Fastener 1: SecondSite™ Adjustable Hub: Catheter Clamp 1: Catheter Trimmer 1: Safety Scalpel: #11 1: Patient ID Card 1: Chart Sticker 1: Patient Information Booklet 1: Flow Rate Information Card with Injection Log 1: Checklist/CLIP Sheet 1: Paper Tape Measure 1: Dressing: STATLOCK®¹ Catheter Stabilization Device |
|--|---|

¹A registered trademark of C. R. Bard, Inc.
All components are CE 0086 unless otherwise noted.



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Fluid path components are non-pyrogenic.

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
CE 0086

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com