



Arterial Catheterization Kit

REF NS-04100-S

20^{Ga.} | 12^{cm} catheter length | .025^{inch dia.} spring-wire guide

Contents:

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- 1: Catheter: 20 Ga. x 1-1/4" (3.2 cm) Jelco^{®1} Catheter over 22 Ga. Introducer Needle
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
- 2: Syringe: 3 mL Luer-Lock
- 1: Drape: 18" x 26" with 3" fenestration
- 1: Prep Sponge Swab
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)

¹ A registered trademark of Medex, Inc.



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Fluid path components are non-pyrogenic.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

DuPont™ Tyvek® packaging

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IS-04100-102D (2/17)



Arterial Catheterization Kit

REF NS-04100-S

20^{Ga.} | 12^{cm} catheter length | .025^{inch dia.} spring-wire guide