

REF NS-04100-S

ARROW



# Arterial Catheterization Kit

20<sup>Ga.</sup> | 12<sup>cm</sup> catheter length | .025<sup>inch dia.</sup> spring-wire guide

**Contents:**

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- 1: Catheter: 20 Ga. x 1-1/4" (3.2 cm) Jelco<sup>®1</sup> Catheter over 22 Ga. Introducer Needle
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
- 2: Syringe: 3 mL Luer-Lock
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Prep Sponge Swab
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)

<sup>1</sup>A registered trademark of Medex, Inc.

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**

**Warning: Cancer and Reproductive Harm** [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Fluid path components are non-pyrogenic.



STERILE EO



Not made with natural rubber latex



REF NS-04100-S

ARROW



# Arterial Catheterization Kit

20<sup>Ga.</sup> | 12<sup>cm</sup> catheter length | .025<sup>inch dia.</sup> spring-wire guide

Manufacturer:  
Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex**<sup>®</sup>  
[www.teleflex.com](http://www.teleflex.com)

IS-04100-105B (2021-03)