



REF ES-04660

24 Ga. 1 Lumen 12 cm catheter length .018 inch dia. spring-wire guide

# One-Lumen CVC

**Contents:**

- 1: One-Lumen Catheter: 24 Ga. (0.8 mm OD) x 12 cm
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Syringe: 3 mL Luer-Slip
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener

All components are CE 2797 unless otherwise noted.

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU) California Prop. 65**

**⚠ WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.**

Fluid path components are non-pyrogenic.



CAS # 7440-48-4



| Lumen           | Priming Volume* (mL) | Gravity Flow Rate† (mL/hr) |
|-----------------|----------------------|----------------------------|
| Distal (24 Ga.) | 0.1                  | 436                        |

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer:

**EC REP** Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
2797 Co. Westmeath, Ireland

**CE** Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL055707 R01 (2023-09)

LOT



24 Ga. 1 Lumen 12 cm catheter length .018 inch dia. spring-wire guide

REF ES-04660

# One-Lumen CVC