

DO NOT BEND

ARROW®



REF SAC-00520-PBX

20^{Ga.}

5^{cm}
catheter
length

Arterial Catheterization Set

Contents:

- 1: Catheter: 20 Ga. x 5 cm Radiopaque PEBA
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW

All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.17 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

en bg cs da nl et fi fr de el hu
 it lv lt no pl pt ro ru sk sl es
 sv tr

LOT



EU Authorized Representative and Importer:

EC REP Teleflex Medical IDA Business and Technology Park Dublin Road, Athlone Co. Westmeath, Ireland

CE 2797 Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

Teleflex® www.teleflex.com

LBL058185 R00 (2022-01)



REF SAC-00520-PBX

20^{Ga.}

5^{cm}
catheter
length

Arterial Catheterization Set

ARROW®

DO NOT BEND