

DO NOT BEND

ARROW®



REF SAC-00822-PBX  
(IPNXXXXXX)

22 Ga. 8 cm catheter length

Arterial Catheterization Set

- Contents:
- 1: Catheter: 22 Ga. x 8 cm Radiopaque PEBA
  - 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
  - 1: Introducer Needle: 22 Ga. x 1-1/2" (4 cm) XTW
- All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.17 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.  
www.teleflex.com/IFU  
California Prop. 65  
⚠ **WARNING:** Cancer and Reproductive Harm.  
www.P65Warnings.ca.gov  
**Not made with natural rubber latex.**  
Fluid path components are non-pyrogenic.




en cs da nl et fi fr de el hu it  
lt no pl pt ro ru sk sl es sv tr

Product of xxx  
Packaged in xxx


LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number



EU Authorized Representative and Importer:

EC REP Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
2797 Co. Westmeath, Ireland

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
www.teleflex.com

LBL058215 R01 (2024-01)



REF SAC-00822-PBX

22 Ga. 8 cm catheter length

Arterial Catheterization Set

ARROW®

DO NOT BEND