

DO NOT BEND



REF SAC-01220

(IPNXXXXXX)

20<sup>Ga.</sup> 12<sup>cm catheter length</sup>

# Arterial Catheterization Set

**Contents:**

- 1: Catheter: 20 Ga. x 12 cm Radiopaque Polyurethane
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 20 Ga. x 2-3/4" (7 cm) XTW

All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.20 mL Priming volume is approximate.

**Rx only Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

[www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.**  
Fluid path components are non-pyrogenic.



- en bg cs da nl et fi fr de el hu
- it lv lt no pl pt ro ru sk sl es
- sv tr

Product of xxx  
Packaged in xxx

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

EU Authorized Representative and Importer:

EC REP Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex**  
[www.teleflex.com](http://www.teleflex.com)

LBL058224 R01 (2024-04)



REF SAC-01220

20<sup>Ga.</sup> 12<sup>cm catheter length</sup>

# Arterial Catheterization Set



DO NOT BEND