

DO NOT BEND

ARROW®



REF SAC-01220-PBX
(IPNXXXXXX)

20^{Ga.} 12^{cm catheter length}

Arterial Catheterization Set

- Contents:**
- 1: Catheter: 20 Ga. x 12 cm Radiopaque PEBA
 - 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
 - 1: Introducer Needle: 20 Ga. x 2-3/4" (7 cm) XTW
- All components are CE 2797 unless otherwise noted.**

Catheter Priming Volume: 0.20 mL **Priming volume is approximate.**

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU
California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



en bg cs da nl et fi fr de el hu
it lv lt no pl pt ro ru sk sl es
sv tr

Product of xxx
Packaged in xxx

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number



EU Authorized
Representative and
Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL058227 R01 (2024-01)



REF SAC-01220-PBX

20^{Ga.} 12^{cm catheter length}

Arterial Catheterization Set

ARROW®

DO NOT BEND